Extreme Comfort™ Particulate Respirator

All products on this page are

MADE
IN U.S.A.

3230+ SURGICAL N95 RESPIRATOR

Up to 50% lower breathing resistance compared to many molded N95 designs.

NIOSH N95 provides at least 95% filtration efficiency against certain non-oil based particles.

3230+ FDA Cleared Surgical N95 Respirator TC-84A-9317

This respirator has been approved as a NIOSH N95 filtering facepiece respirator, for use in healthcare settings, as a Surgical N95 Respirator conforming to recognized standards for biocompatibility, flammability, and fluid resistance.

EXTREME COMFORT™


- Extremely low breathing resistance of 4.1 mm H₂O inhalation and 5.1 mm exhalation - you hardly notice you're wearing an N95 respirator
- Latex-free synthetic rubber head straps for comfortable tension profile
- Staple-free smoothness
- No eyewear fogging
- Soft inner lining

HIGHEST QUALITY

- FDA Cleared, NIOSH Approved Surgical N95 Respirator
- Fluid Resistant – 120 mm HG
- Meets CDC guidelines for TB exposure control
- Individually packaged for hygienic handling

TRUSTED

- GERSON® has been manufacturing NIOSH Approved and FDA Cleared respirators for more than 35 years.
- ISO 9001: 2015 certified manufacturing of respirator protective devices
- Made in Middleboro, MA

Individually
Wrapped
Packages

3230+ SURGICAL N95 RESPIRATOR

DESCRIPTION	QUANTITY #	MODEL #	PART #
Surgical N95 Respirator	50/box, 4 boxes/case	3230+	083230+

GERSON®

Molded vs. Foldable Comparison Chart

MOLDED

Typical Molded N95 (non-valved)

FOLDABLE

3230 Series N95 (non-valved)

		✓ NIOSH Approved Various	✓ NIOSH Approved TC-84A-9285
Protection	#3230 NIOSH Approved N95 Respirator	✓ NIOSH Approved Various	✓ NIOSH Approved TC-84A-9285
	#3230+ Approved for Medical Use	X / ✓	TC-84A-9317
Construction	Durable outer shell	✓	✓
	Low profile, provides optimal field of vision	✓	✓
	Staple free construction	X / ✓	✓
	Nosewire enclosed	X	✓
Comfort	Comfort Level	Moderate	Extreme Comfort
	Up to 50% lower breathing resistance compared to many molded N95 designs.	X	✓
	Extremely reduces eyewear fogging	X	✓
Convenience	Folded, pouch style design for easy carrying and compact storage	X	✓
	Individually packaged to help maintain cleanliness and reduce the risk of contamination.	X	✓
	Easy to carry	X	✓
	Vending-Ready	X	✓

Molded vs. Foldable Testimonials

Testimonials

MOLDED	VS.	FOLDABLE
<i>"Molded N95s are what I used to wear. They are effective BUT..."</i>		<i>"I tried Gerson's Foldable N95! This is the mask I will be wearing in 2021 because..."</i>
<i>"I get face fatigue"</i>		<i>"it is extremely comfortable"</i>
<i>"the straps can be tight"</i>		<i>"the straps are great!"</i>
<i>"the straps are a problem"</i>		<i>"the straps have a nice give and take balance"</i>
<i>"it is difficult to wear all day"</i>		<i>"I can wear it all day - no problem!"</i>
<i>"the dome easily crushes - can't take it anywhere"</i>		<i>"before I need to use it, the folded mask fits in my pocket"</i>

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
1. Unwrap respirator and bend nosepiece into a curve to open mask.
2. Separate the pre-scored elastic headband into two individual bands.
3. Hold the headbands with your index fingers and thumbs.
4. Cup the open mask **securely** under your chin.
5. Pull both headbands over your head.
6. Rest the top headband on the crown of your head and bring the bottom one down under your ears.
7. Press the nosepiece around the contour of your nose and cheekbones. Adjust the edges of the mask to make sure there are no gaps.
8. Perform a Fit Check by gently cupping both hands around the mask, then exhale sharply to see if you sense or detect any leaks.

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WARNING: A respirator does not eliminate the risk of contracting any disease or infection. Occupational use of this respirator requires that you be fit tested by a qualitative or quantitative fit test and that your employer comply with OSHA 29 CFR 1910.134, which includes medical evaluation, training and fit testing.

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100% polypropylene
100% polypropylène
100% polipropileno

WARNINGS

Read and understand warnings
on insert before using.

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© 2021 Louis M. Gerson Co., Inc. All Rights Reserved