

N95 Extreme Comfort™ Particulate Respirator

All products on this page are

MADE
IN U.S.A.



3230 - N95 RESPIRATOR

Up to 50% lower breathing resistance compared to many molded N95 designs.

The Emergency Use Authorization (EUA) issued by the Food and Drug Administration (FDA) on March 2, 2020 permits use of this NIOSH certified respirator in medical settings to mitigate transmission of COVID-19*.

- Durable outer layer and super-soft liner - Made to withstand hardworking environments while providing all day comfort.
- Folded, pouch style design for easy carrying and compact storage.
- Individually packaged to help maintain cleanliness and reduce the risk of contamination.
- Synthetic rubber head straps designed for a comfortable and secure tension profile and extended life. Not made with natural rubber latex.
- Staple-free construction for comfort and safety.
- Adjustable enclosed nosepiece helps reduce potential fogging of eyewear.
- Meets CDC Guidelines for TB Exposure Control.
- NIOSH Approved TC-84-9285.
- Fluid Resistant – 120mm HG.



Individually
Wrapped
Packages

3230 - N95 RESPIRATOR

DESCRIPTION	QUANTITY #	MODEL #	PART #
N95 Respirator	50/box, 4 boxes/case	3230	083230