

GERSON®

1740 - N95 Particulate Respirator



1740 with Exhalation Valve

Larger Comfort-Maxx™ valve opens more easily to expel hot, humid breath and keep inside of the respirator cooler and more comfortable

N95 Particulate Respirators provide at least 95% filtration efficiency against solid and liquid aerosols that do not contain oil. TC-84A-0161

- Soft foam nose cushion
- Contoured formable aluminum nose-clip
- Smooth latex-free polymer headstraps won't catch on hair. No pre-stretching needed
- Safer staple-free construction
- More durable, double-shell construction for long life and better value

PACKAGING							
ITEM NUMBER	NIOSH APPROVAL	MASTER PACKAGING	CASE WEIGHT	CASE DIMENSIONS	CASE CUBE	INNER BARCODE	CASE BARCODE
1740	N95	10/box, 100/Case	6lbs	27" x 11" x 7.5"	1.29 cf	0-44585-81740-0	5-00-44585-91740-2

WARNINGS: Do not use in or around sandblasting. This product does not supply oxygen. Read and Understand Instructions and Warnings before using this product.



GERSON®

16 Commerce Blvd.
Middleboro, MA 02346
Tel 508.947.4000 Fax 508.947.5442
sales@gersonco.com • www.GersonCo.com
© 2021 Louis M. Gerson Co., Inc. All Rights Reserved
Gerson-1740 FeatureSheet MAR21

